

Labycare-Band saw for sectioning body and limbs (SKSS971)



Band saw for sectioning body and limbs

Description : Band saw for sectioning body and limbs

Specifications :

- Max. Throat in mm 200.
 - Table size (L x W) in mm 600 x 600. Blade speed in mtr/min 20-100
 - Blade size (L x W x T) in mm 3505 x 27x0.9.
 - Saw motor capacity in HP 3
 - Coolant motor capacity in HP 0.16.
 - Overall size (LxWx H) in mm 900 x 700 x 2100. Approx Weight in kg 400.
-

