

Labycare-Storage tank to hold cadavers, static/movable, durable tank with input and output facility with lid (SKSS972)



Storage tank to hold cadavers, static/movable, durable tank with input and output facility with lid

Description :

Storage tank to hold cadavers, static/movable, durable tank with input and output facility with lid

Specifications :

Size : 7'5" x 4' 4" x 3

- Technical Specifications For Storage Tank to Hold 10 Cadavers Made of Concrete with Copper Lid with Tilting
- Storage tank to hold 10 cadavers, made of concrete with copper lid Type 304 stain less-steel construction
Built-in cadaver submersion storage system with an.integral cadaver platform that can be raised when needed for study and lowered for storage.
- Plastic Tanks for storing soft and Dissected Part
- Technical Specifications For Plastic Tank For Storing Soft And Dissected Parts

- Should be &best quality.
- Should be a medical grade product.

labycare

www.labycare.com | Email: labycare@gmail.com