

Labycare-Meat cutting machine (SKSS973)



Meat cutting machine

Description :

Meat cutting machine

Specifications :

- Useful for preparing specimen f or big size in Anatomy and Meat Departments. Fitted with large moving table and extension table operated on four ball-bearing rollers. Used extensively in the meat packing, and wholesale fish industry, for handling sword fish and large halbutes etc.
- | | Specifications | | | | (Approx. | |
|--------|----------------|----|-----------|--------|----------|-------|
| Size) | | | | | | |
| | Size | of | cutting | table | 785 | x 585 |
| mm | | | | | | |
| | Total | | table | travel | | |
| 1245mm | | | | | | |
| | | | Extension | table | 455 | x |
| 760mm | | | | | | |
| | Size | of | wheel | | | 455 |
| mm | | | | | | |

- Height 1700 mm
- The table is made of thick S.S. sheet with special heavy axles for easy and firm movement.

Supplied complete with one blade, starter, cord and plug, works on 220 V, single phase, 50Hz , AC supply.

labycare

www.labycare.com | Email: labycare@gmail.com