

## Labycare-Paraffin embedding Bath (SKSS974)



Temperature range from ambient temperature to 95°C. Thermostatic control with an accuracy

of  $\pm 10$ C. Double walled, inside stainless steel and outside mild steel sheet painted in epoxy powder coating.

Top of the bath and concentric rings of stainless steel. To work on 220/230 volts A.C. with standard double

walled water bath but provided with cups & sleeves for glass tubes. 6 cups & 8 sleeves.

### **Description :**

Paraffin embedding Bath

### **Specifications :**

Temperature range from ambient temperature to 95°C. Thermostatic control with an accuracy

of  $\pm 10$ C. Double walled, inside stainless steel and outside mild steel sheet painted in epoxy powder coating.

walled water bath but provided with cups & sleeves for glass tubes. 6 cups & 8 sleeves.

