

Labycare-Ice Lined Refrigerator (SKSS231)



Ice Lined Refrigerator manufacturer by **Labycare** in India.

Description : **Labycare Ice lined refrigerator** is a kind of refrigerator, comprises of high quality ice packs, which help in maintaining 0°C to 8°C temperature without power supply. Such refrigerator is used to carry vaccines, medicines and blood bags.

Specifications :

- Refrigerant: Environment friendly CFC-free
- Temperature Controller: Automatic/Microprocessor
- Temperature display: externally readable digital temperature display
- Power Supply (V/Hz)- 220-240v, 50Hz
- Minimum holdover time at 43°C: 32 hrs
- Eco friendly 134a refrigerant is used. CFC free puff insulation is used.
- Compressor is robust and designed to be used in tropical condition.
- The equipment comes with Solid lid and handle with lock for security.

labycare

www.labycare.com | Email: labycare@gmail.com