

Labycare-Face Shield Disposable 400 Micron (SKSS-CP12)



Face Shield Disposable 400 Micron Manufactured by LABYCARE India

Description :

Face shield Disposable Standard is part of Personal Protective Equipment (PPE) which aims to protect the entire face from hazards such as flying objects, chemical splashes in laboratories and hospitals Face Shield Disposable Standard ensures safety during treatment, diagnosis and quarantine of patients with the novel Coronavirus disease COVID-19

Face Shield Material : 400 Micron thick clear strong robust sheet.

- * Soft Comfortable Grey Foam.
- * Good Quality 1 inch elastic headband of black colour.
- * Both side protective film lamination with removable film to protect with scratches.
- * Odour less & materials used to bring comfort to user.
- * Low fogging face shield.
- * Enough space to wear your spectacles and mask inside face shield.
- * Distance of Shield from Face : 15-20 MM.
- * Coverage : Eyes, Nose, Mouth.

Specifications :

Face shield Disposable Standard is part of Personal Protective Equipment (PPE) which aims to protect the entire face from hazards such as flying objects, chemical splashes in laboratories and

hospitals Face Shield Disposable Standard ensures safety during treatment, diagnosis and quarantine of patients with the novel Coronavirus disease COVID-19

Face Shield Material : 400 Micron thick clear strong robust sheet.

- * Soft Comfortable Grey Foam.
- * Good Quality 1 inch elastic headband of black colour.
- * Both side protective film lamination with removable film to protect with scratches.
- * Odour less & materials used to bring comfort to user.
- * Low fogging face shield.
- * Enough space to wear your spectacles and mask inside face shield.
- * Distance of Shield from Face : 15-20 MM.
- * Coverage : Eyes, Nose, Mouth.

labycare

www.labycare.com | Email: labycare@gmail.com