

Labycare-Spine Board 2 Folded (SKSS1110)



Spine Board 2 Folded Manufacturer By LABYCARE

Description :

This spine board with safe belts is made of PE material with no discharge contaminator, and firmness to wear.

- * This plastic spine board also can be used as a floatable device.
- * This spine board is four-folded, easy for operation and storage.
- * The structure is allowed X-rays because it's completely translucent
- * It is compatible with most head immobilization devices and strap mechanisms.
- * It is mainly used for hospitals, sports, ambulance and outdoors to carry patients and wounded person.

Technical Parameter :

- * Product Size (L x W x H) :- 185 x 50 x 4cm
- * Folding Size :- 94 x 50 x 7 cm
- * Packing Size (1 pc) :- 98 x 10 x 56 cm
- * Load Bearing :- ?159 kg
- * N.W (1 pc) :- 9.2 kg
- * GW :- 10.5kg

We are Exporters and Suppliers of Spine Board 2 Folded at wholesale competitive prices from India.

Specifications :

This spine board with safe belts is made of PE material with no discharge contaminator, and firmness to wear.

- * This plastic spine board also can be used as a floatable device.
- * This spine board is four-folded, easy for operation and storage.
- * The structure is allowed X-rays because it's completely translucent
- * It is compatible with most head immobilization devices and strap mechanisms.
- * It is mainly used for hospitals, sports, ambulance and outdoors to carry patients and wounded person.

Technical Parameter :

- * Product Size (L x W x H) :- 185 x 50 x 4cm
- * Folding Size :- 94 x 50 x 7 cm
- * Packing Size (1 pc) :- 98 x 10 x 56 cm
- * Load Bearing :- 159 kg
- * N.W (1 pc) :- 9.2 kg
- * GW :- 10.5kg

We are Exporters and Suppliers of Spine Board 2 Folded at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com