

Labycare-Oxygen Cylinder Aluminum (A) 5 Ltr. (SKSS-AP059)



Oxygen Cylinder Aluminum (A) 5 Ltr. Manufacturer By Labycare

Description :

A complete oxygen medical kit that is 45 % lighter than traditional steel cylinders.

Refillable and reusable.

With its cylinder made from aluminum alloy, OxyGo is completely rust free unlike steel.

The simple finger operated valve and regulator that can be used by the patient.

Oxygen quality within the cylinder maintained due to its internal aluminum finish and special internal cleaning process.

Working Duration: 6hr 35 min*

Filled Weight of kit: 6 kg

* Duration specified is @2 liter per minute

O2 Capacity : 789 Ltr.

Water Capacity : 5 Ltr.

We are Exporters and Suppliers of Oxygen Cylinder Aluminum (A) 5 Ltr. at wholesale competitive prices from India.

Specifications :

A complete oxygen medical kit that is 45 % lighter than traditional steel cylinders.

Refillable and reusable.

With its cylinder made from aluminum alloy, OxyGo is completely rust free unlike steel.

The simple finger operated valve and regulator that can be used by the patient.

Oxygen quality within the cylinder maintained due to its internal aluminum finish and special internal cleaning process.

Working Duration: 6hr 35 min*

Filled Weight of kit: 6 kg

* Duration specified is @2 liter per minute

O2 Capacity : 789 Ltr.

Water Capacity : 5 Ltr.

We are Exporters and Suppliers of Oxygen Cylinder Aluminum (A) 5 Ltr. at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com