

Labycare-Anaesthesia Machine Silver (SKSS-AP075)



Anaesthesia Machine Silver Manufacturer By Labycare

Description :

Trolley is made of tubular section

Capacious drawer for keeping anaesthesia kit /accessories.

Extended platform at rear for two 10 lits. (Water Cap.) cylinders

Large roller bearings, antistatic castors, for easy mobility and brakes.

Top tray is provided at eye level for keeping monitoring equipments

Hook for breathing circuit.

Wipe clean stainless steel cup.

Rota meter block for oxygen & nitrous oxide with fine adjustment controls & colour coated knob and luminous back plate. Oxygen flow meter calibrated from 100cc/min to 10 ltrs/min

Nitrous oxide flow meter form 200cc/min to 12 litrs/min down

One gold man type vaporizer for halothane.

Non return cum blow off valve in circuit ensures that pressure does not build-up beyond 200 cm.

Water column & minimizes risk of back flow of gases.

60 psl outlet for driving ventilator

Change over system open circuit or close circuit

Emergency Oxygen flush valve at table level to give 35-60 ltrs/min for emergencies

Cylinder yokes : Pin indexed yokes two each of oxygen & nitrous oxide

63 mm diameter colour coded gauges fitted at convenient angle for better visibility & clarity

High efficiency regulators for reliability two each for O₂ & N₂O

One Bain circuit

Safety Features

Emergency oxygen flush

Alarm to provide audible indication of O₂ supply pressure failure

Nitrous oxide shut off in the event of O₂ Failure imension

We are Manufacturer, Exporters and Suppliers of Anaesthesia Machine Silver at wholesale competitive prices from India.

Specifications :

Trolley is made of tubular section

Capacious drawer for keeping anaesthesia kit /accessories.

Extended platform at rear for two 10 lits. (Water Cap.) cylinders

Large roller bearings, antistatic castors, for easy mobility and brakes.

Top tray is provided at eye level for keeping monitoring equipments

Hook for breathing circuit.

Wipe clean stainless steel cup.

Rota meter block for oxygen & nitrous oxide with fine adjustment controls & colour coated knob and luminous back plate. Oxygen flow meter calibrated from 100cc/min to 10 ltrs/min

Nitrous oxide flow meter form 200cc/min to 12 litres/min down

One gold man type vaporizer for halothane.

Non return cum blow off valve in circuit ensures that pressure does not build-up beyond 200 cm.

Water column & minimizes risk of back flow of gases.

60 psl outlet for driving ventilator

Change over system open circuit or close circuit

Emergency Oxygen flush valve at table level to give 35-60 ltrs/min for emergencies

Cylinder yokes : Pin indexed yokes two each of oxygen & nitrous oxide

63 mm diameter colour coded gauges fitted at convenient angle for better visibility & clarity

High efficiency regulators for reliability two each for O₂ & N₂O

One Bain circuit

Safety Features

Emergency oxygen flush

Alarm to provide audible indication of O₂ supply pressure failure

Nitrous oxide shut off in the event of O₂ Failure in-spiration

We are Manufacturer, Exporters and Suppliers of Anaesthesia Machine Silver at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com