

Labycare-Anesthesia Workstation Diamond (SKSS-AP083)



Anesthesia Workstation Diamond Manufacturer By LABYCARE

Description :

Anesthesia Workstation with its ergonomic design, integrates a concise and smooth structure and high mobility.

The microcomputer-controlled Ventilator adopts exported sensors and other components. Its scientific and proper design as well as rigorous production process are what demonstrate the inner quality of the Product.

Mechanical ventilation is made secure and visible by the folded respiratory bellows and the security of air supply is guaranteed by pneumatic transmission and measurement system. With multiple modern technologies applied, the Anesthesia Workstation is both outstanding in performance and versatile in functions.

The Ventilator adopts microcomputer control with high-definition 10.4" TFT color LCD display, integrating the functions of respiratory control and displaying the monitoring parameters. It is versatile in function and excellent in performance.

The high-precision Vaporizer is designed with a stable concentration output, automatic compensation for pressure, temperature and flow rate. And it is available to choose two kinds of gas from Halothane, Enflurane, Isoflurane, and Sevoflurane, ensuring a successful anesthesia operation.

High-definition five tube flowmeter (oxygen, nitrous oxide, and air) are designed with stoppers for oxygen and nitrous oxide to control the output oxygen concentration.

Reliable pneumatic transmission system.

Part of Anesthesia Workstation

Main Body : High-strength engineering plastic rack, light, beautiful and corrosion resistant

Scope of application : Adult and Child

Gas source : O₂, N₂O, Air: 0.28~0.6MPa

Flow meter : O₂: 0.05~1.0L/min 1.1~10L/min

N₂O: 0.05~1.0L/min 1.1~10L/min

Air : 0.05~10L/min

O₂, N₂O linkage and N₂O Stopper

When using nitrous oxide, oxygen concentration>25%;

When the oxygen pressure<0.2Kpa, the flow of nitrous oxide would be cut off.

Specifications :

Anesthesia Workstation with its ergonomic design, integrates a concise and smooth structure and high mobility.

The microcomputer-controlled Ventilator adopts exported sensors and other components. Its scientific and proper design as well as rigorous production process are what demonstrate the inner quality of the Product.

Mechanical ventilation is made secure and visible by the folded respiratory bellows and the security of air supply is guaranteed by pneumatic transmission and measurement system. With multiple modern technologies applied, the Anesthesia Workstation is both outstanding in performance and versatile in functions.

The Ventilator adopts microcomputer control with high-definition 10.4" TFT color LCD display, integrating the functions of respiratory control and displaying the monitoring parameters. It is versatile in function and excellent in performance.

The high-precision Vaporizer is designed with a stable concentration output, automatic compensation for pressure, temperature and flow rate. And it is available to choose two kinds of gas from Halothane, Enflurane, Isoflurane, and Sevoflurane, ensuring a successful anesthesia operation.

High-definition five tube flowmeter (oxygen, nitrous oxide, and air) are designed with stoppers for oxygen and nitrous oxide to control the output oxygen concentration.

Reliable pneumatic transmission system.

Part of Anesthesia Workstation

Main Body : High-strength engineering plastic rack, light, beautiful and corrosion resistant

Scope of application : Adult and Child

Gas source : O₂, N₂O, Air: 0.28~0.6MPa

Flow meter : O₂: 0.05~1.0L/min 1.1~10L/min

N₂O: 0.05~1.0L/min 1.1~10L/min

Air : 0.05~10L/min

O₂, N₂O linkage and N₂O Stopper

When using nitrous oxide, oxygen concentration>25%;

When the oxygen pressure<0.2Kpa, the flow of nitrous oxide would be cut off.

labycare

www.labycare.com | Email: labycare@gmail.com