

Labycare-Laryngeal Mask Silicone Single Use (SKSS-AP085)



Laryngeal Mask Silicone Single Use Manufacturer By LABYCARE

Description :

A laryngeal mask airway is also known as laryngeal mask

It keeps a patient's airway open during anaesthesia or unconsciousness. It is a type of supraglottic airway.

A laryngeal mask is composed of an airway tube that connects to an elliptical mask with a cuff which is inserted through the patient's mouth, down the windpipe, and once deployed forms an airtight seal on top the glottis (unlike tracheal tubes which pass through the glottis) allowing a secure airway.

They are most commonly used by anaesthetists to channel oxygen or anaesthesia gas to a patient's lungs during surgery and in the pre-hospital setting (for instance by paramedics and emergency medical technicians) for unconscious patients

They are for single use & latex free.

We are Manufacturer, Exporters and Suppliers of Laryngeal Mask Silicone Single Use at wholesale competitive prices from India

Specifications :

A laryngeal mask airway is also known as laryngeal mask

It keeps a patient's airway open during anaesthesia or unconsciousness. It is a type of supraglottic airway.

A laryngeal mask is composed of an airway tube that connects to an elliptical mask with a cuff which is inserted through the patient's mouth, down the windpipe, and once deployed forms an airtight seal

on top the glottis (unlike tracheal tubes which pass through the glottis) allowing a secure airway.

They are most commonly used by anaesthetists to channel oxygen or anaesthesia gas to a patient's lungs during surgery and in the pre-hospital setting (for instance by paramedics and emergency medical technicians) for unconscious patients

They are for single use & latex free.

We are Manufacturer, Exporters and Suppliers of Laryngeal Mask Silicone Single Use at wholesale competitive prices from India

labycare

www.labycare.com | Email: labycare@gmail.com