

Labycare-Antistatic Face Mask (SKSS-AP088)



Antistatic Face Mask Manufacturer By LABYCARE

Description :

The Antistatic Face Mask is designed to prevent the buildup of static electricity during medical or industrial procedures. It provides effective respiratory protection while minimizing electrostatic discharge risks. Made from specialized materials, it ensures safety in sensitive environments such as operation theatres, laboratories, and electronic or explosive-prone workplaces.

Specifications :

Antistatic Face Mask are suitable for use in respiratory care of anaesthesia and resuscitator

It is made of antistatic black rubber & come with S. S. Hook & Ring.

It has Air filled cushion to provide safety and comfortable use

These masks have excellent grip with complete air tightness

Electrical resistance is the special feature of Antistatic Mask

Masks have adequate hardness and elasticity

There is no gas leakage

We are Manufacturer, Exporters and Suppliers of Antistatic Face Mask at wholesale competitive prices from India.

