

Labycare-CPR Mask Pocket Type (SKSS-AP093)



CPR Mask Pocket Type Manufacturer By LABYCARE

Description :

The CPR Mask pocket type is a compact, foldable device designed for safe rescue breathing during emergencies. It features a one-way valve to prevent cross-contamination and ensures effective ventilation. Easy to carry in first aid kits, it is widely used by responders, medical staff, and trained individuals.

Specifications :

CPR is known as Cardiopulmonary resuscitation.

It is an emergency procedure that combines chest compression often with artificial ventilation in an effort to manually preserve intact brain function until further measures are taken to restore spontaneous blood circulation and breathing in a person who is in cardiac arrest. It is indicated in those who are unresponsive with no breathing or abnormal breathing, for example, agonal respirations.

Folding cushion mask for patient comfort.

Oxygen inlet for the therapy of oxygen supply in an emergency.

15mm OD connector to use with the manual resuscitator.

Disposable ,detachable one-way valve with 3M filter to prevent cross contamination.

We are Manufacturer, Exporters and Suppliers of CPR Mask Adult Silicone Mask at wholesale competitive prices from India.
