

## Labycare-Oxygen Mask (SKSS-AP094)



Oxygen Mask Manufacturer & Supplier By LABYCARE

### **Description :**

An oxygen mask provides a method to transfer breathing oxygen gas from a storage tank to the lungs.

It is made of clear, non-toxic PVC, medium concentration with Adjustable nose clip

It is Non-Autoclavable.

It comes with oxygen tube.

We are Manufacturer, Exporters and Suppliers of Oxygen Mask at wholesale competitive prices from India.

### **Specifications :**

An oxygen mask provides a method to transfer breathing oxygen gas from a storage tank to the lungs.

It is made of clear, non-toxic PVC, medium concentration with Adjustable nose clip

It is Non-Autoclavable.

It comes with oxygen tube.

We are Manufacturer, Exporters and Suppliers of Oxygen Mask at wholesale competitive prices from India.

