

Labycare-Venturi Mask with 6 Diluters (SKSS-AP098)



Venturi Mask with 6 Diluters Manufacturer & Supplier BY LABYCARE

Description :

Venturi Mask With 6 Diluters is known as Venturi Mask With Six Adopter

The venturi mask is also known as an air-entrainment mask which is used to deliver oxygen concentration to patients on controlled oxygen therapy. This mask is a replacement for intermittent oxygen treatment.

Venturi masks are considered high-flow oxygen therapy devices. This is because venturi masks are able to provide total inspiratory flow at a specified FiO₂ to patients therapy. The kits usually include multiple jets, which are usually color-coded, in order to set the desired FiO₂

It is made of clear, non-toxic PVC, Adjustable nose clip & Non-Autoclavable.

Six color coded diluters: blue 24%, yellow 28%, white 31%, green 35%, pink 40%, orange 50%.

Humidification cup: 22mm-15cm

Flex tube: 20m-22m

Features:

Clear

Soft mask

Straight connector

We are Manufacturer, Exporters and Suppliers of Venturi Mask with 6 Diluters at wholesale competitive prices from India.

Specifications :

Venturi Mask With 6 Diluters is known as Venturi Mask With Six Adopter

The venturi mask is also known as an air-entrainment mask which is used to deliver oxygen concentration to patients on controlled oxygen therapy. This mask is a replacement for intermittent oxygen treatment.

Venturi masks are considered high-flow oxygen therapy devices. This is because venturi masks are able to provide total inspiratory flow at a specified FiO₂ to patients therapy. The kits usually include multiple jets, which are usually color-coded, in order to set the desired FiO₂

It is made of clear, non-toxic PVC, Adjustable nose clip & Non-Autoclavable.

Six color coded diluters: blue 24%, yellow 28%, white 31%, green 35%, pink 40%, orange 50%.

Humidification cup: 22mm-15cm

Flex tube: 20m-22m

Features:

Clear

Soft mask

Straight connector

We are Manufacturer, Exporters and Suppliers of Venturi Mask with 6 Diluters at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com