

## Labycare-Laryngoscope Fibre Optic Miller 4 Blades SS Rexine Pouch (SKSS-AP112)



Laryngoscope Fibre Optic Miller 4 Blades SS Rexine Pouch Manufacturer By LABYCARE

### **Description :**

Fiber Optic Laryngoscope Set has stainless steel construction, a removable light pipe and smooth surface design.

Fiber Optic Laryngoscope Set includes Four Miller blades and small handle

Miller blades are also known as Straight Blades.

Premium quality materials are used in manufacturing of the Laryngoscope Fiber Optic handle and blades.

Single-piece stainless steel construction ensures blade integrity

Removable light pipe and smooth surface design

The free fibre-optic repolishing service extends product life

These sets are custom made combination of blades.

The Blade size available are as below :

Infant Blades : 00, 0

Child Blades : 1, 2

Adult Blades : 3, 4

We are Exporters and Suppliers of Laryngoscope Fibre Optic Miller 4 Blades SS Rexine Pouch at wholesale competitive prices from India.

**Specifications :**

Fiber Optic Laryngoscope Set has stainless steel construction, a removable light pipe and smooth surface design.

Fiber Optic Laryngoscope Set includes Four Miller blades and small handle

Miller blades are also known as Straight Blades.

Premium quality materials are used in manufacturing of the Laryngoscope Fiber Optic handle and blades.

Single-piece stainless steel construction ensures blade integrity

Removable light pipe and smooth surface design

The free fibre-optic repolishing service extends product life

These sets are custom made combination of blades.

The Blade size available are as below :

Infant Blades : 00, 0

Child Blades :1,2

Adult Blades : 3,4

We are Exporters and Suppliers of Laryngoscope Fibre Optic Miller 4 Blades SS Rexine Pouch at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)