

Labycare-Blood Bank Refrigerated Centrifuge 6 and 12 Bags (SKSS879)



Refrigerated centrifuge for 6 and 12 bags ensures precise blood component separation with digital control, efficient cooling, and reliable performance.

Description :

The Blood Bank Refrigerated Centrifuge (6 and 12 Bags) is designed for efficient separation of blood components under controlled temperature conditions. It ensures reliable performance, maintaining sample integrity through precise speed and temperature control. Ideal for blood banks and laboratories, this centrifuge offers quiet operation, digital controls, and durable construction for consistent, high-quality processing.

Specifications :

Specification	Details
Model Capacity	6 Bags / 12 Bags
Type	Refrigerated Centrifuge
Maximum Speed	Up to 5000 RPM
Maximum RCF	Up to 5000 × g
Temperature Range	-20°C to +40°C
Temperature Accuracy	±1°C
Control System	Microprocessor-based digital control
Display	LCD / LED digital display
Timer	Up to 99 minutes
Refrigeration System	CFC-free compressor with rapid cooling
Rotor Type	Swing-out rotor for blood bags
Material	Stainless steel inner chamber
Safety Features	Imbalance detection, lid lock, auto shut-off
Power Supply	220–240V AC, 50Hz
Application	Blood banks, hospitals, research laboratories

labycare

www.labycare.com | Email: labycare@gmail.com