

Labycare-Oxygen Concentrator Dual (SKSS-155)



Oxygen Concentrator Dual

Description :

PSA (Pressure Swing Adsorption) technology

Adjustable flowrate

Large LCD display for: accumulating timing :

Present timing: presetting timing;

Digital: Pressure Temperature

Purity digital; SPO₂, digital

Maintenance reminding

Nebulizer

Flowrate: 0-5L/min, Purity: 93%±3%

Outlet Pressure(Mpa): 0.04-0.07 Sound level (dB): <50

Power: AC 220V±10%, 50/60 ±1Hz;

Power: <480W

Timer: accumulating timing; present timing

Alarm: Power failure, High & Low pressure, High temperature

Low Purity alarm: When oxygen purity is above 92%, it will give green light; When the purity is between 85%,

it will give yellow light; when purity is below 73%, it will give red light;

Maintenance alarm: Alarm for maintenance will be heard after 3000 hours of operation.

We are Manufacturer, Exporters and Suppliers of Oxygen Concentrator Dual Flow at wholesale competitive prices from India.

Specifications :

PSA (Pressure Swing Adsorption) technology

Adjustable flowrate

Large LCD display for: accumulating timing :

Present timing: presetting timing;

Digital: Pressure Temperature

Purity digital; SPO₂, digital

Maintenance reminding

Nebulizer

Flowrate: 0-5L/min, Purity: 93%±3%

Outlet Pressure(Mpa): 0.04-0.07 Sound level (dB): <50

Power: AC 220V±10%, 50/60 ±1Hz;

Power: <480W

Timer: accumulating timing; present timing

Alarm: Power failure, High & Low pressure, High temperature

Low Purity alarm: When oxygen purity is above 92%, it will give green light; When the purity is between 85%,

it will give yellow light; when purity is below 73%, it will give red light;

Maintenance alarm: Alarm for maintenance will be heard after 3000 hours of operation.

We are Manufacturer, Exporters and Suppliers of Oxygen Concentrator Dual Flow at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com