

Labycare-Mox Regulator (SKSS-AP165)



Mox Regulator Manufacturer by LABYCARE

Description :

Mox Regulator is used for A Type Cylinder.

Mox Regulator is made of Brass Forging body with high quality diaphragm

The Regulator has a factory preset fixed static outlet pressure and fitted safety valve.

The regulator automatically maintain flow rate.

It has high quality Pressure Gauge, Rubber Diaphragm, Pressure Gauge which indicates Cylinder pressure and Pressure Gauge which indicates Regulator Outlet Pressure.

We are Manufacturer, Exporters and Suppliers of Mox Regulator at wholesale competitive prices from India.

Specifications :

Mox Regulator is used for A Type Cylinder.

Mox Regulator is made of Brass Forging body with high quality diaphragm

The Regulator has a factory preset fixed static outlet pressure and fitted safety valve.

The regulator automatically maintain flow rate.

It has high quality Pressure Gauge, Rubber Diaphragm, Pressure Gauge which indicates Cylinder pressure and Pressure Gauge which indicates Regulator Outlet Pressure.

We are Manufacturer, Exporters and Suppliers of Mox Regulator at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com