

Labycare-Plasma Extractor Expressor (SKSS880)



SKSS is the largest manufacturer and exporter of Plasma Extractor in India.

The Plasma extractor is designed for separation of blood components from the blood bag.

Plasma Extractor features

- (1) Available in Stainless steel and powder coated mild steel cabinet withstand corrosion, wear & tear on continuous usage.
- (2) The Compression plate is designed to exert uniform pressure on the blood bag.
- (3) SKSS plasma extractor comes with attractive front panel bezel.

Description :

SKSS is the largest manufacturer and exporter of Plasma Extractor in India.

The Plasma extractor is designed for separation of blood components from the blood bag.

Plasma Extractor features

- (1) Available in Stainless steel and powder coated mild steel cabinet withstand corrosion, wear & tear on continuous usage.
- (2) The Compression plate is designed to exert uniform pressure on the blood bag.
- (3) SKSS plasma extractor comes with attractive front panel bezel.

Specifications :

We offer micro-controller based plasma extractors, designed for automatic separation of blood

components from the blood bag. The Plasma Extractor is designed...

Usage : For Laboratory

Capacity : 1 Kg/hr

Material : Iron & Bakelite

Weight : 1000 Kilograms (kg)

labycare

www.labycare.com | Email: labycare@gmail.com