

Labycare-Bougie Straight Catheter (SKSS-AP182)



Bougie Straight Catheter Manufacturer By LABYCARE

Description :

Bougie Straight Catheter has low point of pressure tip which reduces risk of accidental trauma

It has highly visible colour for improved positioning in poor light

The extra tip graduations allow for improved laryngeal inlet entry and exit guidance

The length available in Adult size is 600mm and 800mm length, each with angled (coudé) tip

These Catheters contain barium tip for improved X-Ray guidance (adult sizes only)

Bougie Angled Catheter are latex Free

We are Manufacturer, Exporters and Suppliers of Bougie Straight Catheter at wholesale competitive prices from India.

Specifications :

Bougie Straight Catheter has low point of pressure tip which reduces risk of accidental trauma

It has highly visible colour for improved positioning in poor light

The extra tip graduations allow for improved laryngeal inlet entry and exit guidance

The length available in Adult size is 600mm and 800mm length, each with angled (coudé) tip

These Catheters contain barium tip for improved X-Ray guidance (adult sizes only)

Bougie Angled Catheter are latex Free

We are Manufacturer, Exporters and Suppliers of Bougie Straight Catheter at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com