

Labycare-Semi Electric OT Table General (SKSS-OT206)



Semi Electric OT Table General Manufacturer by LABYCARE

Description :

A Semi Electric OT Table is designed for general surgical procedures, combining electric motorized movements with manual controls for reliable and cost-effective operation.

Specifications :

Operation: Electric height adjustment; manual Trendelenburg, lateral tilt & back/leg sections

Tabletop: Radiolucent, X-ray compatible

Height Range: Approx. 700–1000 mm

Trendelenburg / Reverse: $\pm 25^\circ$

Lateral Tilt: $\pm 15-20^\circ$

Load Capacity: 180–250 kg

Control: Handset + manual levers

Base: Stainless steel base with central locking castors

Power Supply: 220–240 V AC

Accessories: Head rest, arm boards, leg supports, anesthesia screen

labycare

www.labycare.com | Email: labycare@gmail.com