

Labycare-Electric C-Arm Table with Top Slide (SKSS-OT209)



Electric C-Arm Table with Top Slide

Description :

An Electric C-Arm Table with top slide is designed for orthopedic, trauma, and vascular procedures, allowing smooth longitudinal movement of the tabletop for enhanced C-arm imaging without repositioning the patient.

Specifications :

Operation: Fully electric (height, Trendelenburg, lateral tilt, top slide)

Tabletop: Radiolucent, X-ray compatible sliding top

Top Slide Travel: Approx. 300–350 mm

C-Arm Access: Unobstructed head-to-toe and lateral access

Height Range: 700–1000 mm (approx.)

Trendelenburg / Reverse: $\pm 25^\circ$

Lateral Tilt: ± 15 – 20°

Load Capacity: 200–250 kg

Control: Handset / remote control

Base: Stainless steel base with central locking system

Power Supply: 220–240 V AC, 50/60 Hz

Standard Accessories: Arm boards, leg supports, headrest, anesthesia screen

labycare

www.labycare.com | Email: labycare@gmail.com