

Labycare-Ophthalmic OT Table Manual (SKSS-OT211)



Ophthalmic OT Table Manual

Description :

A Manual Ophthalmic OT Table is specially designed for eye surgeries, providing smooth manual adjustments, stable positioning, and patient comfort during delicate ophthalmic procedures.

Specifications :

- **Operation:** Fully manual (height, backrest & Trendelenburg as applicable)
 - **Tabletop:** Cushioned, antistatic, radiolucent top
 - **Head Section:** Special ophthalmic headrest with fine adjustment
 - **Height Range:** Approx. 750–1000 mm
 - **Backrest Adjustment:** Manual crank / lever operated
 - **Load Capacity:** 150–200 kg
 - **Base & Frame:** Stainless steel construction
 - **Mobility:** Antistatic castors with brakes
 - **Applications:** Cataract, retinal, glaucoma & other eye surgeries
-

labycare

www.labycare.com | Email: labycare@gmail.com