

Labycare-Ophthalmoscope Welch Allyn 12860 (SKSS-OI217)



Ophthalmoscope Welch Allyn 12860

Description :

Welch Allyn PocketScope Ophthalmoscope Set, Part No. 12860 (Single)

Lightweight, compact diagnostic instruments deliver high performance in a slightly smaller package.

Halogen light for true tissue colour and consistent, long-lasting illumination

12 aperture/filter combinations for greater versatility: micro, small, and large spot sizes, cobalt blue filter for corneal exams, fixation target, and slit aperture, plus red-free filter and unfiltered halogen

48 focusing lenses for better resolution

2,5V PocketScope battery handle is durable, lightweight and small enough to carry in your pocket. Features an adjustable rheostat that controls light intensity and a compact pocket clip. Works with two AA batteries.

We are Exporters and Suppliers of Ophthalmoscope Welch Allyn 12860 at wholesale competitive prices from India.

Specifications :

Welch Allyn PocketScope Ophthalmoscope Set, Part No. 12860 (Single)

Lightweight, compact diagnostic instruments deliver high performance in a slightly smaller package.

Halogen light for true tissue colour and consistent, long-lasting illumination

12 aperture/filter combinations for greater versatility: micro, small, and large spot sizes, cobalt blue filter for corneal exams, fixation target, and slit aperture, plus red-free filter and unfiltered halogen

48 focusing lenses for better resolution

2,5V PocketScope battery handle is durable, lightweight and small enough to carry in your pocket. Features an adjustable rheostat that controls light intensity and a compact pocket clip. Works with two AA batteries.

We are Exporters and Suppliers of Ophthalmoscope Welch Allyn 12860 at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com