

Labycare-Haemocytometer Marienfeld (SKSS-PLP229)



Haemocytometer Marienfeld

Description :

The haemocytometer is a counting-chamber device used for counting blood cells. Haemocytometer set comes in a plastic case consisting of :

1 counting chamber without clamps with double net ruling.

Two cover glasses,

One blood diluting pipette Thoma for red blood cells.

One blood diluting pipette Thoma for white blood cells

Two silicone tubings of approx. 16 cm length

One red plastic mouth piece

One white plastic mouthpiece

We are Exporters and Suppliers of of Haemocytometer Marienfeld at wholesale competitive prices from India.

Specifications :

The haemocytometer is a counting-chamber device used for counting blood cells. Haemocytometer set comes in a plastic case consisting of :

1 counting chamber without clamps with double net ruling.

Two cover glasses,

One blood diluting pipette Thoma for red blood cells.

One blood diluting pipette Thoma for white blood cells

Two silicone tubings of approx. 16 cm length

One red plastic mouth piece

One white plastic mouthpiece

We are Exporters and Suppliers of of Haemocytometer Marienfeld at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com