

Labycare-Haemocytometer Standard (SKSS-PLP230)



Haemocytometer Standard Manufacturer & Supplier

Description :

The hemocytometer is a counting-chamber device used for counting blood cells. Haemocytometer set comes in a plastic case consisting of :

1 counting chamber without clamps with double net ruling.

Two cover glasses,

One blood diluting pipette Thoma for red blood cells.

One blood diluting pipette Thoma for white blood cells

Two silicone tubings of approx. 16 cm length

One red plastic mouth piece

One white plastic mouthpiece

We are Manufacturer, Exporters and Suppliers of Haemocytometer Standard at wholesale competitive prices from India.

Specifications :

The hemocytometer is a counting-chamber device used for counting blood cells. Haemocytometer set comes in a plastic case consisting of :

1 counting chamber without clamps with double net ruling.

Two cover glasses,

One blood diluting pipette Thoma for red blood cells.

One blood diluting pipette Thoma for white blood cells

Two silicone tubings of approx. 16 cm length

One red plastic mouth piece

One white plastic mouthpiece

We are Manufacturer, Exporters and Suppliers of Haemocytometer Standard at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com