


# Labycare-Counting Chamber Neubauer Marienfeld (SKSS-PLP231)



Counting Chamber Neubauer Marienfeld Manufacturer & Supplier

## **Description :**

Counting Chamber is made of Thick Glass.

Counting Chamber is used for blood cell counting, the squares placed at the corners are used for white cell counting. Since their concentration is lower than red blood cells a larger area is required to perform the cell count. Therefore, the central square is made of 400 small squares.

We are Exporters and Suppliers of Counting Chamber Neubauer Marienfeld at wholesale competitive prices from India.

## **Specifications :**

Counting Chamber is made of Thick Glass.

Counting Chamber is used for blood cell counting, the squares placed at the corners are used for white cell counting. Since their concentration is lower than red blood cells a larger area is required to perform the cell count. Therefore, the central square is made of 400 small squares.

We are Exporters and Suppliers of Counting Chamber Neubauer Marienfeld at wholesale competitive prices from India.

---

labycare

www.labycare.com | Email: labycare@gmail.com