


Labycare-Counting Chamber Neubauers Standard (SKSS-PLP232)



Counting Chamber Neubauers Standard

Description :

Counting Chamber is made of Thick Glass.

Counting Chamber is used for blood cell counting, the squares placed at the corners are used for white cell counting. Since their concentration is lower than red blood cells a larger area is required to perform the cell count. Therefore, the central square is made of 400 small squares.

We are Manufacturer, Exporters and Suppliers of Counting Chamber Neubauers Standard at wholesale competitive prices from India.

Specifications :

Counting Chamber is made of Thick Glass.

Counting Chamber is used for blood cell counting, the squares placed at the corners are used for white cell counting. Since their concentration is lower than red blood cells a larger area is required to perform the cell count. Therefore, the central square is made of 400 small squares.

We are Manufacturer, Exporters and Suppliers of Counting Chamber Neubauers Standard at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com