


## Labycare-R.B.C. Pipette Standard (SKSS-PLP234)



R.B.C. Pipette Standard Manufacturer By LABYCARE

### **Description :**

RBC pipette is diluting pipette.

RBC pipette is known as Red Blood Cell Pipette.

The blood along with the RBC diluting fluid is drawn in the pipette and are mixed well through the red bead present inside the bulb of the pipette.

The size of the bulb is larger in RBC pipette when compared to the size WBC pipette.

For Red Blood Corpuscles (Erythrocytes, Proportion Of Dilution 1 : 100)

White Back

Blue Graduation

Tolerance:  $\pm 3 \%$

Divided Into 10 Parts

We are Manufacturer, Exporters and Suppliers of R.B.C. Pipette Standard at wholesale competitive prices from India.

### **Specifications :**

RBC pipette is diluting pipette.

RBC pipette is known as Red Blood Cell Pipette.

The blood along with the RBC diluting fluid is drawn in the pipette and are mixed well through the red

bead present inside the bulb of the pipette.

The size of the bulb is larger in RBC pipette when compared to the size WBC pipette.

For Red Blood Corpuscles (Erythrocytes, Proportion Of Dilution 1 : 100)

White Back

Blue Graduation

Tolerance:  $\pm 3\%$

Divided Into 10 Parts

We are Manufacturer, Exporters and Suppliers of R.B.C. Pipette Standard at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)