

Labycare-W.B.C. Pipette Standard (SKSS-PLP236)



W.B.C. Pipette Standard Manufacturer & Supplier

Description :

WBC pipette is diluting pipette.

WBC pipette is known as White Blood Cell Pipette.

WBC diluting fluid is used for performing the WBC (Leucocyte) count.

Draw EDTA anticoagulated blood to 0.5 mark in the capillary end of WBC pipette and Draw diluting fluid up to 11 mark.

For White Blood Corpuscles (Leucocytes, Proportion Of Dilution 1 : 10)

White Back

Blue Graduation

Tolerance: $\pm 3 \%$

Divided Into 10 Parts

We are Manufacturer, Exporters and Suppliers of W.B.C. Pipette Standard at wholesale competitive prices from India.

Specifications :

WBC pipette is diluting pipette.

WBC pipette is known as White Blood Cell Pipette.

WBC diluting fluid is used for performing the WBC (Leucocyte) count.

Draw EDTA anticoagulated blood to 0.5 mark in the capillary end of WBC pipette and Draw diluting fluid up to 11 mark.

For White Blood Corpuscles (Leucocytes, Proportion Of Dilution 1 : 10)

White Back

Blue Graduation

Tolerance: $\pm 3\%$

Divided Into 10 Parts

We are Manufacturer, Exporters and Suppliers of W.B.C. Pipette Standard at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com