

## Labycare-ESR Tube Glass Marienfeld (SKSS-PLP244)



ESR Tube Glass Marienfeld Manufacturer By LABYCARE

### **Description :**

ESR Tube is made of Borosilicate Glass.

ESR Tube is also known as Westergren sedimentation Tube.

ESR tubes measure the rate at which red blood cells settle out of plasma.

Westergren sedimentation tubes are used for the standard procedure, by drawing blood in a liquid sodium citrate vacuum tube

Westergren tube is open at both ends.

It is 30 cm in length and 2.5 mm in diameter.

The lower 20 cm are marked with 0 at the top and 200 at the bottom.

We are Exporters and Suppliers of ESR Tube Glass Marienfeld at wholesale competitive prices from India.

### **Specifications :**

ESR Tube is made of Borosilicate Glass.

ESR Tube is also known as Westergren sedimentation Tube.

ESR tubes measure the rate at which red blood cells settle out of plasma.

Westergren sedimentation tubes are used for the standard procedure, by drawing blood in a liquid sodium citrate vacuum tube

Westregren tube is open at both ends.

It is 30 cm in length and 2.5 mm in diameter.

The lower 20 cm are marked with 0 at the top and 200 at the bottom.

We are Exporters and Suppliers of ESR Tube Glass Marienfeld at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)