

Labycare-Stainless Steel Stretcher on Trolley (SKSSS378B)



Stainless Steel Stretcher on Trolley manufacturer by **Labycare** in India.

Description :

Labycare Stainless Steel Stretcher on Trolley is often equipped with variable height frames, wheels, tracks, or skids. Stretchers are primarily used in acute out-of-hospital care situations by emergency medical services (EMS), military, and search and rescue personnel.

Specifications :

- Frame work made of 1.25" x 16G vertical and 1" x 18G horizontal
- CRC tubes
- Trolley mounted on 15 cms. dia castors, 2 with brakes
- Removable stretchers top made of M.S. Sheet with S.S. handle
- Overall size: 213 L x 56 W x 81 H cms
- Finish: Pre-treated and Epoxy Powder Coated
- Optional
- Mattress
- Side railings M.S./S.S

labycare

www.labycare.com | Email: labycare@gmail.com