

Labycare-Manual ICU Bed 3 Function (SKSS801K)



The Manual Patient Bed is designed to provide safety, comfort, and reliable support for patients during treatment and recovery. It features smooth manual crank operation for easy positioning of backrest, footrest, and Trendelenburg functions. Built with a durable epoxy-coated structure and patient-friendly components, the bed ensures stability, hygiene, and ease of use in hospitals, nursing homes, and clinical environments.

Description :

Labycare Manual ICU Bed

Backed by the team of expert professionals we are engaged in manufacturing, supplying and exporting ICU Bed. Offered beds are known for their robust construction, comfort level, adjustable height and durability. These beds are made in strict compliance with the international standards. Once our ICU Beds are manufactured, these are tested on various parameters under the supervision of our expert team. Customers can avail these beds from us at market leading prices.

Specifications :

Technical Specifications

- **Bed Type:** Manual, 3-function (3 crank system)
- **Functions:** Backrest, footrest, Trendelenburg adjustment
- **Operation:** Manual crank mechanism
- **Frame Material:** Epoxy-coated mild steel
- **Bed Platform:** 3-section mattress support
- **Head & Foot Board:** ABS, detachable with locking system
- **Side Rails:** Collapsible ABS with angle indicator
- **Mobility:** 125 mm noiseless castors with individual brakes
- **Safety Features:** Corner bumper protection
- **Accessories:** Adjustable IV pole, drainage hooks
- **Overall Dimensions:** 2160 × 960 × 500 mm
- **Backrest Adjustment:** 0° – 75°
- **Footrest Adjustment:** 0° – 40°

- **Trendelenburg:** 0° – 12°
 - **Safe Working Load:** 140 – 240 kg
 - **Standards & Certification:** CE Mark, MDD 93/42/EEC, ISO 13485
-

labycare

www.labycare.com | Email: labycare@gmail.com