

Labycare-Patient Transfer Trolley (SKSSHS378A)



Patient Transfer Trolley manufacturer by **Labycare** in India. We have established ourselves as a renowned supplier, manufacturer and exporter of superior quality of Patient transfer trolley.

Description :

Labycare Patient Trolley is used for carrying injured patients from one place to another. Offered trolley's frame is made up of stainless steel tubes that are very strong. Over this frame, a stretcher is placed which can be removed from the frame, if required. Our manufactured trolley is made in strict compliance with the international standards. Moreover, Patient Transfer Trolley is also equipped with the removable side safety railings.

Specifications :

- Frame work made of 1.25" x 16G vertical and 1" x 18G horizontal
- CRC tubes
- Trolley mounted on 15 cms. dia castors, 2 with brakes
- Removable stretchers top made of M.S. Sheet with S.S. handle
- Overall size: 213 L x 56 W x 81 H cms
- Finish: Pre-treated and Epoxy Powder Coated

Optional

- Mattress
- Side railing M.S./S.S.

labycare

www.labycare.com | Email: labycare@gmail.com