


Labycare-Test Tube Stand Round 12 Holes (SKSS-PLP272)



Test Tube Stand Round 12 Holes

Description :

Test Tube Stand Round by Plastic Labware India is a handy little space saver with circular twelve place Polypropylene .

Test Tube Stand has four 25 mm dia holes & eight 19 mm dia holes in its top plate.

The base has matching hemispherical wells & 16 vertical pins for drying tubes surrounded by a trough to catch drippings.

This autoclavable rack can be easily disassembled for cleaning purpose.

We are Manufacturer, Exporters and Suppliers of Test Tube Stand Round 12 Holes at wholesale competitive prices from India.

Specifications :

Test Tube Stand Round by Plastic Labware India is a handy little space saver with circular twelve place Polypropylene .

Test Tube Stand has four 25 mm dia holes & eight 19 mm dia holes in its top plate.

The base has matching hemispherical wells & 16 vertical pins for drying tubes surrounded by a trough to catch drippings.

This autoclavable rack can be easily disassembled for cleaning purpose.

We are Manufacturer, Exporters and Suppliers of Test Tube Stand Round 12 Holes at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com