

Labycare-Head Immobilizer (SKSSHI215)



Head Immobilizer Supplier by **Labycare** in India.

Description :

Labycare Head Immobilizer fixes a patient's head in place to prevent further injury after suffering trauma. The head immobilizer is made of high-density, plastic material and contains large ear holes for monitoring the patient's ear canal.

Specifications :

- This Head Immobilizer is a new generation product of our company.
- It is made from 3 pcs separate part, light weight, wasy to carry and use, made of NBR Rubber Foam.
- It can be used in the winter and summer, afford the impact of temperature difference.
- There are two holes which are used to observe the bleeding or drainage on both sides of the head, hole size is 6.5 cm in diameter.
- *The material of product does not contain metal.
- The patient can take product during examination of X-Ray, CT, nuclear magnetic resonance imaging (MRI).

labycare

www.labycare.com | Email: labycare@gmail.com