

## Labycare-Clot Activator Tube Non Vacuum (SKSS-PLP290)



Clot Activator Tube Non Vacuum MANUFACTURER BY LABYCARE

### **Description :**

Clot Activator Tube is Coated with micronised Silica particles which causes blood to clot rapidly.

Clot Activator Tube provides enough and non polluted serum specimen for clinical test.

Cap Colour : Red

Max. Vacuum Degree Tolerance : +10 %

Invert the Tube for 6-8 times gently and centrifuge for 10 minutes after 15-30 minutes coagulation.

It is used in Blood Collection and storage for Biochemistry, Immunology and Serology tests.

We are Manufacturer, Exporters and Suppliers of Clot Activator Tube Non Vacuum at wholesale competitive prices from India.

### **Specifications :**

Clot Activator Tube is Coated with micronised Silica particles which causes blood to clot rapidly.

Clot Activator Tube provides enough and non polluted serum specimen for clinical test.

Cap Colour : Red

Max. Vacuum Degree Tolerance : +10 %

Invert the Tube for 6-8 times gently and centrifuge for 10 minutes after 15-30 minutes coagulation.

It is used in Blood Collection and storage for Biochemistry, Immunology and Serology tests.

We are Manufacturer, Exporters and Suppliers of Clot Activator Tube Non Vacuum at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)