

## Labycare-Gel & Clot Activator Tube Vacuum (SKSS-PLP291)



Gel & Clot Activator Tube Vacuum

### **Description :**

Gel & Clot Activator Tube Vacuum

Gel & Clot Activator Tube is Coated with micronised Silica particles which causes blood to clot rapidly and the barrier gel can effectively separate the serum from fibrin and cells, while preventing substance exchange between blood cell and serum.

Owing to its clot activator, it is better for keeping the original character of the blood and solving the problem of hemolysis.

Cap Colour : Golden Yellow

Max. Vacuum Degree Tolerance : +10 %

Invert the Tube for 6-8 times gently and centrifuge for 10 minutes after 3000-3500rpm after 15-30 minutes coagulation.

It is used in Blood Collection and storage for Biochemistry, Immunology and Serology tests.

We are Manufacturer, Exporters and Suppliers of Gel & Clot Activator Tube Vacuum at wholesale competitive prices from India.

### **Specifications :**

Gel & Clot Activator Tube is Coated with micronised Silica particles which causes blood to clot rapidly and the barrier gel can effectively separate the serum from fibrin and cells, while preventing substance exchange between blood cell and serum.

Owing to its clot activator, it is better for keeping the original character of the blood and solving the

problem of hemolysis.

Cap Colour : Golden Yellow

Max. Vacuum Degree Tolerance : +10 %

Invert the Tube for 6-8 times gently and centrifuge for 10 minutes after 3000-3500rpm after 15-30 minutes coagulation.

It is used in Blood Collection and storage for Biochemistry, Immunology and Serology tests.

We are Manufacturer, Exporters and Suppliers of Gel & Clot Activator Tube Vacuum at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)