

Labycare-Aluminium Scoop Stretcher (SKSSASS216)



Aluminium Scoop Stretcher manufacturer by **Labycare** in India. We are a renowned and trustworthy firm in this domain, engaged in offering an excellent quality of Aluminium Scoop Stretcher.

Description :

Labycare Scoop Stretchers are the ideal stretchers for those who have spinal injuries or injured people. The Scoop reduces the need for the log-roll maneuver, which significantly decreases cervical spine displacement. Two linked, interlocking sections allow operators to bring the two halves under the patient mutually and gently raise them. This stretcher is especially designed to carry patients in hospitals and nursing homes. The offered stretcher is manufactured by our adept professionals utilizing high quality basic materials with the help of innovative techniques at our latest production set up. The provided range of Aluminium Scoop Stretcher is available with us in different sizes.

Specifications :

- Easy to move
 - Lightweight
 - Comfortable
 - High load bearing capacity
 - Aluminium Structure
 - Gross weight Scoop Stretcher: 11 kg
 - Net Weight Of Scoop Stretcher : 9 kg
 - The bearing capacity of Scoop Stretcher : <159 kg
 - Scoop Stretcher size (L x W x H) : 168 x 43 x 7 cm (min) 224 x 45 x 10 cm (max)
 - Dimension of the packing box : 127 x 45 x 10 cm
-

