

Labycare-Contrast Bath (SKSS-PP316)



Contrast Bath Manufacturer by labycare

Description :

Suitable for treatment of arms and legs.

Two stainless steel heavy gauge tanks with inside measurement of 50cm X 38cm X 40cm deep are mounted on 10cm size ball bearing rubber castors.

One tank is fitted with special heating unit and other with heavy duty cooling unit.

Each tank is provided with digital Temp. Controller, pilot lights and water outlets.

Both tanks are well insulated and placed side by side in a heavy tubular steel frame covered with painted steel sheet.

Works on 220V.

We are Exporters and Suppliers of Contrast Bath (Heat & Cold) at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com