

Labycare-Infra Red Lamp 150 Watt Floor Model (SKSS-PP326)



Infra Red Lamp 150 Watt Floor Model Manufacturer by LABYCARE

Description :

Floor stand model lamp for IR radiation treatments.

IR Bulb is fitted in metal casing with wire guard, & automatic cut off timer & an intensity control switch.

Infra-Red Bulb : One, IR Bulb of 150w (Phillip or Osram) is fitted in a metal casing with wire guard.

Intensity Control : Fitted with intensity control switch to control out put.

Timer : Automatic cut off Timer of 60 minutes.

Construction : Bulb is mounted on a wide tubular base (25" x 25") with four ball bearing castors.

Placement : Adjustable arm permits wide range of adjustments in Height and Length for proper placement of heating lamp.

Height is adjustable from 45" to 70".

Horizontal adjustment is upto 18" with 360 degree rotational adjustment.

Finish : Powder coated finish, with chrome plated working parts.

We are Exporters and Suppliers of Infra Red Lamp 150 Watt Floor Model at wholesale competitive prices from India.

Specifications :

Floor stand model lamp for IR radiation treatments.

IR Bulb is fitted in metal casing with wire guard, & automatic cut off timer & an intensity control switch.

Infra-Red Bulb : One, IR Bulb of 150w (Phillip or Osram) is fitted in a metal casing with wire guard.

Intensity Control : Fitted with intensity control switch to control out put.

Timer : Automatic cut off Timer of 60 minutes.

Construction : Bulb is mounted on a wide tubular base (25" x 25") with four ball bearing castors.

Placement : Adjustable arm permits wide range of adjustments in Height and Length for proper placement of heating lamp.

Height is adjustable from 45" to 70".

Horizontal adjustment is upto 18" with 360 degree rotational adjustment.

Finish : Powder coated finish, with chrome plated working parts.

We are Exporters and Suppliers of Infra Red Lamp 150 Watt Floor Model at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com