

Labycare-Infra Red Lamp 450 Watt Floor Model (SKSS-PP327)



Infra Red Lamp 450 Watt Floor Model Manufacturer by LABYCARE

Description :

Floor Stand model Lamp 450watts, for IR Radiation Treatments.

Infra-Red Bulb : Three, IR Bulb of 150w (Phillip or Osram) is fitted in a metal casing with wire guard.

Intensity Control : Fitted with intensity control switch to control out put.

Timer : Automatic cut off Timer of 30 minutes.

Construction : Bulb are mounted on a wide powder coated base with lockable ball bearing casters.

Placement : Adjustable arm permits wide range of adjustments in Height and Length.

We Exporters and Suppliers of Infra Red Lamp 450 Watt Floor Model at wholesale competitive prices from India.

Specifications :

Floor Stand model Lamp 450watts, for IR Radiation Treatments.

Infra-Red Bulb : Three, IR Bulb of 150w (Phillip or Osram) is fitted in a metal casing with wire guard.

Intensity Control : Fitted with intensity control switch to control out put.

Timer : Automatic cut off Timer of 30 minutes.

Construction : Bulb are mounted on a wide powder coated base with lockable ball bearing casters.

Placement : Adjustable arm permits wide range of adjustments in Height and Length.

We Exporters and Suppliers of Infra Red Lamp 450 Watt Floor Model at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com