

Labycare-Infra Red Lamp 750 Watt Floor Model (SKSS-PP328)



Infra Red Lamp 750 Watt Floor Model Manufacturer by LABYCARE

Description :

Nickle chromium resistance wire embedded into a ceramic casting provides uniform heat transmission.

96% thermal efficiency with low element surface temperature and no visible light is produced.

Ceramic Heating Element : 750watts. with NO Visible light and fitted in a straight metal casing with wire guard.

Timer : Automatic cut off Timer of 30 minutes.

Heat Up Time : 3 to 9 minutes to reach 90 degree of maximum temperature.

Intensity Control : Fitted with intensity control switch to control out put.

Placement : Adjustable arm permits wide range of adjustment in Height and Length for proper placement of heating element.

Mounted on wide base with four lockable ball bearing casters.

We are Exporters and Suppliers of Infra Red Lamp 750 Watt Floor Model at wholesale competitive prices from India.

Specifications :

Nickle chromium resistance wire embedded into a ceramic casting provides uniform heat transmission.

96% thermal efficiency with low element surface temperature and no visible light is produced.

Ceramic Heating Element : 750watts. with NO Visible light and fitted in a straight metal casing with wire guard.

Timer : Automatic cut off Timer of 30 minutes.

Heat Up Time : 3 to 9 minutes to reach 90 degree of maximum temperature.

Intensity Control : Fitted with intensity control switch to control out put.

Placement : Adjustable arm permits wide range of adjustment in Height and Length for proper placement of heating element.

Mounted on wide base with four lockable ball bearing casters.

We are Exporters and Suppliers of Infra Red Lamp 750 Watt Floor Model at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com