

Labycare-Shortwave Diathermy 500 Watt (SKSS-PP330)



Shortwave Diathermy 500 Watt

Description :

High power medical diathermy incorporating a high power transmitting oscillator valve and solid state rectifiers; provides R.F. power upto 500 watts for intense deep tissue heating without discomfort at the surface. Housed in an elegant MOBILE CABINET & designed to meet heavy work load in a hospital.

Technical :

Frequency : 27.12MHz

Output : Continuous, Maximum : 500 watts.

Intensity : Adjustable in 5 Steps.

Electrodes : One Pair Condenser pads with felt spacer & cover, Size : 15cm X 25cm

Tuning : Manual, with tuning control knob.

Timer : 30 minutes timer with auto cut.

Cooling : With 10cm dia. instrumental fan.

Voltage : 220-240V

Accessories :

(a) One pair condenser pad cable

(b) One power main cable

(c) Spare Fuses

(d) One shortwave test tube

(e) Two Velcro strap to tie / hold condenser pads with the body

(f) One automatic voltage stabilizer.

We are Exporters and Suppliers of Shortwave Diathermy 500 Watt at wholesale competitive prices from India.

Specifications :

High power medical diathermy incorporating a high power transmitting oscillator valve and solid state rectifiers; provides R.F. power upto 500 watts for intense deep tissue heating without discomfort at the surface. Housed in an elegant MOBILE CABINET & designed to meet heavy work load in a hospital.

Technical :

Frequency : 27.12MHz

Output : Continuous, Maximum : 500 watts.

Intensity : Adjustable in 5 Steps.

Electrodes : One Pair Condenser pads with felt spacer & cover, Size : 15cm X 25cm

Tuning : Manual, with tuning control knob.

Timer : 30 minutes timer with auto cut.

Cooling : With 10cm dia. instrumental fan.

Voltage : 220-240V

Accessories :

(a) One pair condenser pad cable

(b) One power main cable

(c) Spare Fuses

(d) One shortwave test tube

(e) Two Velcro strap to tie / hold condenser pads with the body

(f) One automatic voltage stabilizer.

We are Exporters and Suppliers of Shortwave Diathermy 500 Watt at wholesale competitive prices from India.

