

Labycare-Shortwave Diathermy Solid State 500 Watt (SKSS-PP332)



Shortwave Diathermy Solid State 500 Watt Manufacturer by LABYCARE

Description :

Solid-State, Pre-programmed, Pulsed Shortwave diathermy with solid state rectifiers provides R.F. Power up to 250 watts in Continuous mode and 1000 watts peak power in Pulse mode for intense deep tissue heating without discomfort at the surface.

Designed to meet the heavy work load in a busy hospital.

Features :

Output : Continuous Mode 500 watts, adjustable : Pulse mode 1000 watts adjustable

Timer : 00-30 minutes Adjustable

Pulse Rate : 20,50,100,150,200,250 Hz

Pre-Programmed : 20 Programmed fended

Wave Shape : Mono Phonic

Operation Output Frequency : 27.12MHz

Wave Length : 11 meters

Display : LCD Screen display

Input : 220V AV - 50Hz

Accessories :

(a) One pair condenser pad cable

- (b) One power main cable
- (c) Spare Fuses
- (d) One shortwave test tube
- (e) Two Velcro strap to tie / hold condenser pads with the body
- (f) One automatic voltage stabilizer.

We are Exporters and Suppliers of Shortwave Diathermy Solid State 500 Watt at wholesale competitive prices from India.

Specifications :

Solid-State, Pre-programmed, Pulsed Shortwave diathermy with solid state rectifiers provides R.F. Power up to 250 watts in Continuous mode and 1000 watts peak power in Pulse mode for intense deep tissue heating without discomfort at the surface.

Designed to meet the heavy work load in a busy hospital.

Features :

Output : Continuous Mode 500 watts, adjustable : Pulse mode 1000 watts adjustable

Timer : 00-30 minutes Adjustable

Pulse Rate : 20,50,100,150,200,250 Hz

Pre-Programmed : 20 Programmed fended

Wave Shape : Mono Phonic

Operation Output Frequency : 27.12MHz

Wave Length : 11 meters

Display : LCD Screen display

Input : 220V AV - 50Hz

Accessories :

- (a) One pair condenser pad cable
- (b) One power main cable
- (c) Spare Fuses
- (d) One shortwave test tube
- (e) Two Velcro strap to tie / hold condenser pads with the body
- (f) One automatic voltage stabilizer.

We are Exporters and Suppliers of Shortwave Diathermy Solid State 500 Watt at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com