

Labycare-Muscle Stimulator with Tens Economy (SKSS-PP334)



Muscle Stimulator with Tens Economy MANUFACTURER BY LABYCARE

Description :

The IC Based with Large LCD display.

The combination unit produces TENS & Four different types of stimulating current i.e Galvanic, Interrupted Galvanic, Faradic & Surged Faradic for therapeutic use output is smoothly adjustable from 0-135 volts.

The stimulating current outputs are as under:

GALVANIC : Continuous well filtered D.C upto 135 volts for a maximum of 0-75mA load.

INTERRUPTED GALVANIC : Pulse duration - 0.01 , 0.03 , 0.1 , 0.3 , 1,3 , 10 , 30 , 100 , 300m secs, Frequency 6 ,18 ,30 per minute 1 ,3 ,5 ,10 ,30 ,50 ,100 per sec.

SURGED FARADIC : Faradic Pulse of 0.6 m sec duration, Contraction duration : 2-4 seconds.

FARADIC : Comprises of 80 pulses per second, each pulse of 0.6 ms width.

Accessories :

1 Ball Tip Electrode with Handle ;

2 Pair Electrodes Wires.

We are Exporters and Suppliers of Muscle Stimulator with Tens Economy at wholesale competitive prices from India.

Specifications :

The IC Based with Large LCD display.

The combination unit produces TENS & Four different types of simulating current i.e Galvanic, Interrupted Galvanic, Faradic & Surged Faradic for therapeutic use output is smoothly adjustable from 0-135 volts.

The simulating current outputs are as under:

GALVANIC : Continuous well filtered D.C upto 135 volts for a maximum of 0-75mA load.

INTERRUPTED GALVANIC : Pulse duration - 0.01 , 0.03 , 0.1 , 0.3 , 1,3 , 10 , 30 , 100 , 300m secs, Frequency 6 ,18 ,30 per minute 1 ,3 ,5 ,10 ,30 ,50 ,100 per sec.

SURGED FARADIC : Faradic Pulse of 0.6 m sec duration, Contraction duration : 2-4 seconds.

FARADIC : Comprises of 80 pulses per second, each pulse of 0.6 ms width.

Accessories :

1 Ball Tip Electrode with Handle ;

2 Pair Electrodes Wires.

We are Exporters and Suppliers of Muscle Stimulator with Tens Economy at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com