

Labycare-Ultrasound Therapy Unit Digital Pre Programmed (SKSS-PP339)



Ultrasound Therapy Unit Digital Pre Programmed Manufacturer & Supplier

Description :

Micro Processor Based computerized ultrasound therapy unit gives you a choice of treatment modes continuous or in Three Pulsed wave forms.

The compact digital unit with built in timer and user friendly control panel gives a wide range of treatment situations.

Technical :

Ultrasound Frequency : 1 MHz's $\pm$ 10%

Output Mode : Continuous & Pulsed

Pulsed Frequency : Variable, 16Hz, 48Hz, 99Hz

Intensity : Continuous mode 0 to 2 watts / cm sq. , Pulsed Mode 0 to 3 watts / cm sq.

Program : 24 Preset Programmed Parameters (Intensity, Mode & Time of Treatment).

Treatment Head : One. Effective Area Radiation : 5sq. Centimeters.

Timer : 01 to 60 minutes adjustable, with Buzzer at the end of time.

Display : Two Line digital display of parameters

Power Required : 220V AC, 50Hz.

Accessories :

Ultrasound Gel Bottle.

We are Exporters and Suppliers of Ultrasound Therapy Unit Digital Pre Programmed at wholesale competitive prices from India.

Specifications :

Micro Processor Based computerized ultrasound therapy unit gives you a choice of treatment modes continuous or in Three Pulsed wave forms.

The compact digital unit with built in timer and user friendly control panel gives a wide range of treatment situations.

Technical :

Ultrasound Frequency : 1 MHz's+/-10%

Output Mode : Continuous & Pulsed

Pulsed Frequency : Variable, 16Hz, 48Hz, 99Hz

Intensity : Continuous mode 0 to 2 watts / cm sq. , Pulsed Mode 0 to 3 watts / cm sq.

Program : 24 Preset Programmed Parameters (Intensity, Mode & Time of Treatment).

Treatment Head : One. Effective Area Radiation : 5sq. Centimeters.

Timer : 01 to 60 minutes adjustable, with Buzzer at the end of time.

Display : Two Line digital display of parameters

Power Required : 220V AC, 50Hz.

Accessories :

Ultrasound Gel Bottle.

We are Exporters and Suppliers of Ultrasound Therapy Unit Digital Pre Programmed at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com