

Labycare-Stainless Steel Patient's Revolving Stool (SKSSRS863)



Stainless Steel Patient's Revolving Stool manufacturer by **Labycare** in India.

Description :

Labycare Patient's Revolving Stool is extensively used in hospitals, health care centers, and clinics. It is made of the first-rate alloyed metal that is capable of enduring heavy loads , rust, and shear forces. The stool has got a cushioning seat that gives extreme comfort to the patients. It revolves around its axis smoothly and is usually kept near the patient's bed.

Specifications :

- Four leg base fitted with PVC stumps
- Height adjustment by machine screw from 18" x 27"
- Rest ring support made up of M.S. tube
- S.S. top supported with M.S. sheet top
- Finish: Pre-treated and Epoxy Powder Coated
- Also available in S.S. Frame work

labycare

www.labycare.com | Email: labycare@gmail.com