

Labycare-Ultrasound Therapy Unit Digital Single Head (SKSS-PP340)



Ultrasound Therapy Unit Digital Single Head

Description :

Ultrasound Therapy Unit Digital Single Head

Solid State Ultrasound Therapy Unit gives you a choice of Two treatment modes, Continuous or in Four pulsed wave forms.

The compact unit with built in digital timer & user friendly control panel gives a wide range of treatment situation.

Technical :

Ultrasonic Output Frequency : 1 MHz (+/- 10%)

Ultrasonic Output Modes : Two, Continuous & Pulsed.

Ultrasonic Output Power : 15watts in Continuous mode : 21watts in Pulsed mode.

Ultrasonic Output Intensity in w / cm² : 3.5 cm²

Treatment Head : One, 5cm² large size head.

Pulse Ratio : Four, 1:1 , 1:2 , 1:4 , 1:7 , & 1:10

Timer : 0 to 99 minutes , pre-settable.

Buzzer / Alarm : Audio Buzzer at the end of treatment time .

Display : All Parameters are digitally displayed.

Power : 220V AC / 50Hz.

We are Exporters and Suppliers of Ultrasound Therapy Unit Digital Single Head at wholesale competitive prices from India.

Specifications :

Ultrasound Therapy Unit Digital Single Head

Solid State Ultrasound Therapy Unit gives you a choice of Two treatment modes, Continuous or in Four pulsed wave forms.

The compact unit with built in digital timer & user friendly control panel gives a wide range of treatment situation.

Technical :

Ultrasonic Output Frequency : 1 MHz (+/- 10%)

Ultrasonic Output Modes : Two, Continuous & Pulsed.

Ultrasonic Output Power : 15watts in Continuous mode : 21watts in Pulsed mode.

Ultrasonic Output Intensity in w / cm² :3.5 cm²

Treatment Head : One, 5cm² large size head.

Pulse Ratio : Four, 1:1 , 1:2 , 1:4 , 1:7 , & 1:10

Timer : 0 to 99 minutes , pre-settable.

Buzzer / Alarm : Audio Buzzer at the end of treatment time .

Display : All Parameters are digitally displayed.

Power : 220V AC / 50Hz.

We are Exporters and Suppliers of Ultrasound Therapy Unit Digital Single Head at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com